


GETTING STARTED

ViQtor is easy to use. You wear viQtor under your clothing directly on the upper arm. ViQtor is splash-proof and can be worn while showering. ViQtor works both indoors and outdoors. So, you can wear it anytime.



- A** Personal HELP button
- B** Sensor that detects if you have fallen
- C** Sound
- D** Status LED light, primarily shows battery status: red, orange, green)
- E** Button to check battery status

Step 1: The charger

Remove the charger from the viQtor box. The charging cable can be found at the bottom of the box. Plug this cable with the USB-C connector into the back of the charger. Now it is ready for use.

Place the charger in a visible place with limited direct sunshine, for example, in your lounge.

Make sure the monitoring device is placed on the charger so that the elevation on the back of the device falls into the recess of the charger. When placed correctly, you will hear a beep, see a vibration and the LED light will be purplen.

Step 2: The monitoring device (viQtor)

Make sure viQtor monitoring device is charged:

- Briefly press the small button at the bottom of viQtor **E**. A light will now illuminate to indicate the extent to which the battery is charged. When the light turns green, it means the device is more than 20% charged. If the light illuminates an orange colour, then it is less than 20% charged and you should charge the monitoring device.

Step 3: The armband

The armband consists of two parts:

- A strap that can be tightened using the buckle and adjusted to the right size with Velcro.
- A holder into which viQtor is to be placed has an opening at the top with a cap that you pull over the device.
- When viQtor is in the holder, the sensors on the back of viQtor should be visible through the opening at the back of the holder. Note that the indicated triangle on the back of viQtor goes into the holder first.

Step 4: Wearing on the upper arm

Place the armband on the outer side of the upper arm, slightly above the arm muscle. The armband should be sufficiently tight in order to function properly, but still comfortable, so not too tight. Get into the habit of charging the device for at least 30 minutes to an hour every day. For this purpose, pick a specific moment during the day when you are with someone, such as a family member or a caregiver.

